



JollyPop® Pacifier

The Pacifier Used in Top US Children's Hospitals*

Soft, Comfortable & Lightweight

- 100% silicone - made in USA
- Single-piece construction adheres to AAP guidelines.
- Lowest profile – fits all CPAP, nasal cannulas and side lying, without cutting down.

Cutting down pacifiers violates product safety requirements and presents a risk of injury and infection to caregivers and babies.

NEWBORN For babies without teeth

UNSCENTED

VPJP1301100 Teal
VPJP1308100 Pink
VPJP1310100 Blue
VPJP1303100 Orange
VPJP1312100 Lavender

VANILLA SCENTED

VPJP1302100 Teal
VPJP1309100 Pink
VPJP1311100 Blue
VPJP1323100 Orange
VPJP1313100 Lavender

PREEMIE

UNSCENTED

VPJP1304100 Teal
VPJP1306100 Pink
VPJP1307100 Blue

VANILLA SCENTED

VPJP1305100 Teal

PLUS+ For babies beginning to teethe

UNSCENTED

VPJP1314100 Solid Teal
VPJP1315100 Clear



Over
10 million
used in US
hospitals

JollyPop Hospital Trial Results

98%

of nurses approved
switching to JollyPop**
over 39% wanted
to switch ASAP



**of 922 nurses indicating a preference

Specifications: 100% medical grade silicone. Phthalate-free, latex-free, BPA-free, DEHP-free. Packaged clean & ready to use in individual poly bags. Manufactured in the USA.

* As ranked by U.S. News & World Report

CoverMe™ Disposable Covers For Swings and Rockers

Fits mamaRoo® and many other swings with 3-point harness.

Finally a disposable cover that protects swings and rockers and makes clean-up easy!

- Single patient use
- Medical grade, breathable, waterproof fabric.
- Fits the mamaRoo® classic and many other swings and rockers used in the NICU and special care nursery.

"Made clean-up so easy
after baby got messy."
Knoxville, TN

Item #	Description	Qty.
VPMC160125	Disposable cover	1/bag, 25 per box
VPMC16015	Disposable cover	1/bag, 5 per box

Specifications:

Medical grade, breathable, waterproof fabric. 100% polypropylene cover. Latex-free. Fabric meets AAMI PB70 Level 3 testing. Single-patient use, disposable. For hospital and home use. Do not leave baby unattended in swings or rockers. Manufactured in the USA.

Helps Comply
with Infection
Control



Jellies™ Gel Pillow & Positioning Aid

Safe, affordable, long-lasting gel support.

Thicknesses and sizes to meet all developmental needs.

- 6 x 6 inch Pillow in three thickness
Thin: for micro-preemies & preemies
Medium: for full-term newborns & infants
Thick: for children & adults
- 6 x 9 inch Cushion
Ideal mattress for micro-preemies or pressure relief
- Disposable, baby-soft covers available for each size

Item #	Description	Qty. Per Box
VPJE1504212	6 x 6 x 3/8 inch Thin Gel Pillow	12
VPJE1504312	6 x 6 x 5/8 inch Medium Gel Pillow	12
VPJE1504512	6 x 6 x 3/4 inch Thick Gel Pillow	12
VPJE1504412	6 x 9 x 1/2 inch Gel Cushion	12
VPJE1504812	6 x 9 x 3/4 inch Gel Cushion	12
VPJE1504835	6 x 9 x 3/4 inch Gel Cushion, w/ printed cover	35
VPJE15041130	ULTRAGel 6 x 9 x 1 inch Cushion, w/ printed cover	30
VPJE15047100	Disposable printed covers (for 6 x 6 inch)	100
VPJE15046100	Disposable printed covers (for 6 x 9 inch)	100

Philips Gel-E Donut *	Compares to Jellies:
92025-A	Gel-E Donut Premie Thin
92025-B	Gel-E Donut Small Medium
92025-C	Gel-E Donut Thick
91033-2	Squishon 2 6 x 9 inch
91033-3	Squishon 3 6 x 9 inch

*Gel-E Donut is a registered trademark of Philips Medical.



**Try our new
higher viscosity ULTRA Gel!
Ask for samples.**

Specifications:
Single patient use. Latex free.
Outer Material: soft, polyolefin film.
Contents: Water based non-toxic food grade gel.
Methylparaben anti-bacterial agent.
See SDS for more information.
Manufactured in the USA.

Tucker Sling™ Positioning Aid

Support for babies who need elevation.

- One-piece support system for babies who need to be in elevated positions.
- May be used for babies with GER/acid reflux or respiratory problems.
- Washable.
- Item descriptions indicate size of mattress and weight of baby.

The Tucker Sling is the exact same product I introduced at Children's Medical Ventures. Sandbox Medical became the distributor in January 2016 to ensure these products remain available to customers. Tucker Sling is manufactured by the inventor, Terry Jarret.

*Read McCarty, Co-Founder, Sandbox Medical
Founder, Past CEO Children's Medical Ventures*

Item #	Description	Qty. Per Box
VPPTS16030112	Incubator; Infant <2 lbs. Flannel	12
VPPTS16030212	Incubator; Infant 2-5 lbs. Flannel	12
VPPTS16030312	Bassinet; Infant 3-7 lbs. Flannel	12
VPPTS16030412	Crib; Infant 2-5 lbs. Flannel	12
VPPTS16030512	Crib; Infant 3-7 lbs. Flannel	12
VPPTS16030612	Crib; Infant 6-16 lbs. Flannel	12
VPPTS16030712	Crib; Infant 6-16 lbs. Quilted	12
VPPTS16030812	Crib; Infant 16-24 lbs. Quilted	12
VPPTS16030912	Crib; Infant 24-30 lbs. Quilted	12
VPPTS1603101	Foam wedge 12" x 26" x 24"	1

Specifications: 100 % cotton flannel Sling. Soy-based Foam Wedge.
Manufactured in the USA.



CuddleStar™ Pacifier Stabilizer

Helps keep the pacifier in place.

- 7 x 7 inch star-shaped pacifier holder with removable Preemie or Newborn JollyPop Pacifier.
- Soft non-woven fabric – hospital use only.
- Partially filled with polypropylene beans – safely contained in a separate inner sack.
- Velcro® tab makes it easy to replace JollyPop.
- Single patient use.

Item #	Description	Qty. Per Box
VPCS1502112	CuddleStar with JollyPop Preemie Pacifier	12
VPCS1502212	CuddleStar with JollyPop Newborn Pacifier	12

Specifications:

Phthalate-free, latex-free, BPA-free and DEHP-free. Polyester non-woven outer and inner layer. New, polypropylene beads contained in inner layer. Meets ASTM F963-11 safety code requirements. Manufactured in the USA.



CozyCare Bunting™ Thermal Wrap

Helps prevent heat loss and provides containment

For the delivery room, nursery, transport, feeding or at home.

- Made of soft, non-allergenic, non-toxic and flame retardant materials.
- Layer of DuPont™ Thermolite® helps baby maintain temperature.

The NOVAPLUS CozyCare Bunting is the exact same product I introduced at Children's Medical Ventures. SandBox Medical became the distributor in January 2016 to ensure these products remain available to Vizient members.

Item #	Description	Qty.
VPCB1602125	For babies 6 - 15 lbs.	1/bag, 25 per box
VPCB1602225	For babies 3 - 6 lbs.	1/bag, 25 per box
VPCB1602325	For babies 1 - 3 lbs.	1/bag, 25 per box

Specifications:

100% polyester non-woven fabric. Single-patient use, disposable. Manufactured in the USA.



NP5-010116P

Hospital Ordering Information

Contact us for samples or order direct:

T 1-888-986-0042 orders@chsusainc.com www.chsltd.com

Order from your distributor:

- Cardinal
- Medline
- Seneca
- Claflin
- Owens and Minor
- Shared Services