



JelliesTM

Gel Pillow & Positioning Aid



Meeting All Developmental Needs

6" x 6" Pillows

- Thin: for micropreemies & preemies
- Medium: for fullterm babies
- Thick: for children & adults

6" x 9" Cushions

- Where a larger support is required
- Ideal mattress for micropreemies

New UltraGel!

- Higher viscosity 6" x 9" cushions

Single patient use

Disposable, baby-soft covers available for each size

Made in the USA.



Safe, Affordable, Long-Lasting Gel Support
Give every baby the benefit of gel pillows.

Soft, pliable support for baby's head and body.

- Helps to minimize head molding and skin breakdown.
- Helps to relieve pressure and encourage proper positioning.
- Provides uniform support for the soft cranium of the pre-term infant.
- Preservatives in the nontoxic gel eliminate risk of mold and fungus.
- Lasts 4 to 6 weeks under normal conditions.
- Also makes a great hot or cold pack for nursing moms!

Soft, pliable Jellies help reduce the need for frequent repositioning.



Bright blue color won't get lost in the bed.

For more support and for bigger babies, push the gel to one side, fold under and tape.

Research has shown that gel mattresses are especially effective in reducing vibration during transport.*

*Journal of Perinatology (2015) 35, 110–114; doi:10.1038/jp.2014.172; published online 27 November 2014



With Jellies[™], Every Baby Can Benefit from Gel Support That's Affordable & Long-Lasting.

Gel pillows are an effective way to minimize pressure and provide uniform support for the soft cranium of the pre-term infant. The soft surface and supple support of Jellies Gel Pillows keeps babies more comfortable and minimizes the re-positioning needed to alleviate pressure points.

Easy To Use, Find and Store.

With its bright blue color, Jellies won't get lost in the bed, so you'll save on replacement pillows. You can fold Jellies to adjust the size and thickness. So one size can be used as baby grows. That helps to reduce inventory requirements, ordering and stocking.

Product Specifications

Latex free.

Outer Material: soft polyolefin film.

Contents: Water based non-toxic food grade gel. Methylparaben anti-bacterial agent. See SDS for more information.

Instructions For Use

Single Patient Use Only.

Prewarm gel in incubator or by kneading under warm water (in plastic bag). DO NOT microwave. Cover with a soft cloth or disposable cover. See Jellies complete *Instructions For Use* packaged with product.

Warnings

Jellies are indicated as single patient use to prevent possible cross contamination between patients. Any reuse of Jellies is not warranted by the manufacturer.

Hospital use only. Jellies should only be used in a monitored environment.

DO NOT allow infant to lay face down on pillow. Danger of suffocation. DO NOT X-Ray through pillow. DO NOT microwave.

Jellies comes in thicknesses and sizes to meet all developmental needs.

Item #	Description	Qty Per Box
VPJE1504212	6 x 6 x 3/8 inch Thin Gel Pillow	12
VPJE1504312	6 x 6 x 5/8 inch Medium Gel Pillow	12
VPJE1504512	6 x 6 x 3/4 inch Thick Gel Pillow	12
VPJE1504412	6 x 9 x 1/2 inch Gel Cushion	12
VPJE1504812	6 x 9 x 3/4 inch Gel Cushion	12
VPJE1504835	6 x 9 x 3/4 inch Gel Cushion w/ printed cover	35
VPJE15041130	ULTRAGel 6 x 9 x 1 inch Cushion, w/ printed cover	30
VPJE15047100	Disposable printed covers (for 6 x 6 inch)	100
VPJE15046100	Disposable printed covers (for 6 x 9 inch)	100

Philips Gel-E Donut *	Compares to Jellies:
92025-A Gel-E Donut Premie	Thin
92025-B Gel-E Donut Small	Medium
92025-C Gel-E Donut	Thick
91033-2 Squishon 2	6 x 9 inch
91033-3 Squishon 3	6 x 9 inch

*Gel-E Donut is a registered trademark of Philips Medical.

6 x 6 inch square comes in three thicknesses - to meet the needs of micro-preemies, infants, children & adults.



6 x 9 inch rectangular Jellies – perfect mattress for micro-preemies or hip-roll for bigger babies.



NP5-1504-006n