

Jellies®

Gel Pillow & Positioning Aid



Sizes To Meet All Developmental Needs

- 6 x 6 inch Pillow
5/8 or 3/8 inch thickness
- 6 x 9 inch Cushion
1/2 inch thickness
Ideal mattress for micro-preemies

Disposable,
baby-soft covers
Available for each size

Passes all USP
preservative
effectiveness tests.

Made in the USA.



Safe, Affordable, Long-Lasting Gel Support

Give every baby the benefit of gel pillows.

Jellies lasts your baby's entire length of stay!

Helps minimize head molding and skin breakdown

- Guaranteed mold and bacteria free. Salt based gel will not support the growth of mold, fungus or bacteria.
- Provides uniform support for the soft cranium of the pre-term infant
- Lasts significantly longer than other brands
- Ideal for micro-preemies, preemies and full-term babies
- Also makes a great hot or cold pack for nursing moms!
- Single patient use only



Soft, pliable Jellies help reduce the need for frequent repositioning. For more support and for bigger babies, push the gel to one side, fold under and tape.

Jellies' Salt-Based Gel Is Guaranteed Mold, Fungus and Bacteria-Free.

It's well known that salt prevents bacteria from growing due to its high osmolarity. We chose a salt-based gel to ensure Jellies is safe from the inside out. The gel in Jellies is made of water and sodium chloride (salt) so it will not support the growth of mold, fungus or bacteria.

Jellies nonporous film prevents evaporation to ensure single patient use for entire length of stay.

Jellies passes USP preservative effectiveness tests with results indicating fewer mold, fungus and bacteria than in USP purified drinking water.*

*Independent laboratory tests on file.



Contact us at: orders@sandboxmedical.com

With Jellies®, Every Baby Can Benefit from Gel Support That's Affordable & Long-Lasting.

Gel pillows are an effective way to minimize pressure and provide uniform support for the soft cranium of the pre-term infant. The soft surface and supple support of Jellies Gel Pillows keeps babies more comfortable and minimizes the re-positioning needed to alleviate pressure points. Jellies also can aid in the healing process if sores or ulcers are present. The Jellies nonporous film means they last for your baby's entire length of stay... saving you money while giving all babies the benefits of gel support!

Easy To Use, Find and Store.

With its bright blue color, Jellies won't get lost in the bed, so you'll save money on replacement pillows. And it's nonporous film prevents evaporation to ensure single patient use for entire length of stay.*

Product Specifications

Outer Material: polyethylene film.

Contents: sodium chloride, water, pulverized cellulose, sodium benzoate, food grade gum (guar, xanthan gum).

Product is latex-free.

*Jellies are indicated as single patient use to prevent possible cross contamination between patients. Any reuse of Jellies is not warranted by the manufacturer and should be made only after an appropriate method of sanitizing / disinfecting has been validated by the user.

Jellies comes in two thicknesses and sizes to meet developmental needs of all babies.

Item #	Description	Quantity
JE1504312	5/8 x 6 x 6 inch thick Gel Pillow	12 per box
JE1504212	3/8 x 6 x 6 inch thick Gel Pillow	12 per box
JE1504550	Disposable Covers (6 x 6 inch)	50 per box
JE1504412	1/2 x 6 x 9 inch Gel Cushion	12 per box
JE1504650	Disposable Covers (6 x 9 inch)	50 per box

6 x 6 inch square comes in two thicknesses – to meet the needs of micro-preemies and bigger babies.



6 x 9 inch rectangular Jellies – perfect mattress for micro-preemies or hip-roll for bigger babies.



5-1504-006a