



CozyCare Bunting™

Thermal Wrap



Helps prevent heat loss and provides containment.

For the delivery room, nursery, transport, feeding or at home.

- Made of soft, non-allergenic, non-toxic and flame retardant materials.
- Layer of DuPont™ Thermolite® helps baby maintain temperature.
- Easily adjustable with soft, refastenable tabs.
- Available in three sizes.

The NOVAPLUS CozyCare Bunting is the exact same product I introduced at Children's Medical Ventures. SandBox Medical became the distributor in January 2016 to ensure these products remain available to Vizient members.

CozyCare Bunting continues to be manufactured by ALunt Designs.

Read McCarty, Co-Founder, SandBox Medical
 Founder, Past CEO Children's Medical Ventures

Item #	Description	Qty.
VPCB1602125	For babies 6 - 15 lbs.	1/bag, 25 per box
VPCB1602225	For babies 3 - 6 lbs.	1/bag, 25 per box
VPCB1602325	For babies 1 - 3 lbs.	1/bag, 25 per box

Specifications

100% polyester non-woven fabric. Single-patient use, disposable.
 Manufactured in the USA.



NP5-1602-003c